

EasyCare

User and Testing Manual

A clear guide for installing, operating, demonstrating, and manually testing the EasyCare elderly-assistance app.

Large controls
Designed for simple one- or two-tap use.

Voice or touch
Important actions always have button alternatives.

Important: EasyCare is a prototype. Food orders are not sent to a real restaurant. Calls and messages open in the phone's normal apps so the user can review them before pressing Call or Send. In an immediate or life-threatening emergency, call 112 directly.

Prepared for EasyCare users, caregivers, and app testers

1. What you need

- An Android phone with EasyCare installed.
- Internet access for the best voice-recognition experience.
- Microphone permission if you want to use voice commands.
- Location permission if you want to include your location in an emergency message.

2. Install and open EasyCare

1. Install EasyCare on the Android phone.
2. Find the EasyCare icon on the phone.
3. Tap the icon to open the app.
4. The first-time caregiver setup screen will appear.

3. First-time caregiver setup

EasyCare shows the setup screen only on first launch.

1. Enter or review the user's name.
2. Enter the family member's name and phone number.
3. Enter the emergency contact phone number.
4. Enter the home address.
5. Tap Save and Continue.

The prototype starts with these example values:

- User: Test User
- Family member: Test Family Member
- Family phone: +491234567890

Use test phone numbers when demonstrating the app. The saved details remain available after the app is closed and reopened.

4. Home screen

The home screen shows the user's greeting, current date and time, and the main actions:

- Tap and Speak - use a voice command.
- Food - select a mock meal.
- Call Family - open the family member in the phone dialer.
- Message Family - type or dictate an SMS.
- Medicine - view and complete reminders.
- EMERGENCY HELP - open emergency actions.

The small settings icon is intended for a caregiver.

5. Use voice commands

1. From the home screen, tap Tap and Speak.
2. If prompted, allow microphone access.
3. Wait for Listening....
4. Say one short command clearly.
5. EasyCare displays the detected sentence, speaks a response when voice output is enabled, and opens the matching action.

Example phrases:

What you want	Example phrases
Food	"I want food", "Order food", "I am hungry"
Call family	"Call Dipesh", "Call family", "Call my son"
Message family	"Send message to Dipesh", "Message my son", "Tell my family"
Emergency	"Help me", "Emergency", "I need help"
Medicine	"Medicine", "My medicine", "Medication"

For a call command, EasyCare opens the system dialer directly. It does not place the call automatically. Press the phone app's Call button only when you intend to make the call.

If voice recognition is unavailable or permission is denied, use the large buttons on the home screen. All important functions still work without voice commands.

6. Order food

1. Tap Food, or say a food voice command.
2. Review the four local meal choices.
3. Tap Select on a meal.
4. EasyCare immediately shows Order Request Sent with a Preparing status.
5. Return home when finished.

Available mock meals:

- Gujarati Thali - €10
- Khichdi Kadhi - €8
- Punjabi Thali - €12
- Simple Home Meal - €7

No money is charged and no restaurant receives the request.

7. Call a family member

1. Tap Call Family on the home screen, or open Family and tap Call Dipesh.
2. EasyCare opens the system phone app with the saved number.
3. Review the number.
4. Press Call in the phone app only if you want to place the call.

There is no extra confirmation page inside EasyCare. The phone app is the final review step.

8. Send a family message

1. Tap Message Family.
2. Select Family or Emergency Contact as the recipient.
3. Type a message in the large text box, or tap Speak Message and dictate it.
4. Tap Send Message.
5. Review the prepared message in the system SMS app.
6. Press Send in the SMS app when ready.

If microphone access is denied, typing remains available.

9. Use Emergency Help

Tap the large red EMERGENCY HELP button on the home screen. Choose the family member or emergency contact before preparing an alert.

Call family

Tap Call Dipesh. The system dialer opens with the saved number. Press Call in the dialer to place the call.

Send a quick emergency alert

1. Select the intended recipient.
2. Tap Send Quick Alert.
3. Allow location access if you want to include the current location.
4. EasyCare opens the SMS app with a preset help message.
5. If a location is available, the text includes a Google Maps link.
6. Review the recipient and message, then press Send in the SMS app.

The message can still be prepared if location permission is denied or the device cannot obtain a location.

Speak an emergency message

1. Select the intended recipient.
2. Tap Speak Emergency Message.
3. Allow microphone access if prompted.
4. Speak a short message, such as "I fell down and need help."
5. Check the detected text shown by EasyCare.
6. EasyCare opens the SMS app with the dictated text and an optional location link.
7. Review it and press Send in the SMS app.

Call emergency services

Tap Call Emergency Services. EasyCare opens the system dialer with 112. The app never silently places the call; press Call in the dialer to continue.

10. Mark medicine as taken

1. Tap Medicine.
2. Find the reminder you completed.
3. Tap its large TAKEN button.
4. Confirm that the button changes to DONE ✓.

EasyCare says “Medicine marked as taken” when voice output is enabled. This is only a reminder prototype and does not provide medical advice. Completion state resets when the app process is fully restarted.

11. Caregiver settings

Tap the small settings icon on the home screen. A caregiver can change:

- My Name
- Family Contact
- Emergency Contact
- Voice Assistant ON or OFF
- Text Size: Normal, Large, or Extra Large

Turn off Voice Assistant to stop spoken Text-to-Speech feedback. This does not remove the large navigation buttons.

12. Five-minute demonstration script

Use this sequence for a client or caregiver demonstration:

1. Complete first-launch setup and arrive at the home screen.
2. Tap Tap and Speak and say “I want food.”
3. Select Gujarati Thali and show Order Request Sent and Preparing.
4. Return home, tap Tap and Speak, and say “Call Dipesh.”
5. Show the number in the system dialer, then return without placing a real call.
6. Open EMERGENCY HELP, select a contact, and tap Send Quick Alert.
7. Allow location access and show the prepared SMS and Google Maps link. Return without sending it to a real contact.
8. Open Medicine, tap TAKEN, and show DONE ✓.

9. Open settings, change the text size, and turn voice output off and back on.

13. Manual test checklist

Use test contacts and do not press the final system Call or Send button unless the recipient has agreed to participate.

Test	Steps	Expected result
First launch	Clear app data and open EasyCare	Setup screen appears with example values
Save setup	Change a name, tap Save and Continue, close and reopen	Home screen appears and the changed name is retained
Home content	Open home	Greeting, date, time, voice button, four cards, and emergency button are visible
Food voice command	Say "I want food"	Detected sentence appears and Food opens
Call voice command	Say "Call Dipesh"	System dialer opens with the saved family number
Message voice command	Say "Message my son"	Message screen opens
Emergency voice command	Say "I need help"	Emergency screen opens
Medicine voice command	Say "Medication"	Medicine screen opens
Unknown voice command	Say an unrelated sentence	App explains that it did not understand and remains usable
Microphone denied	Deny microphone permission, then try voice	Clear explanation appears; buttons and typing still work; no crash
Food order	Select every meal one at a time	Success screen shows the selected meal and Preparing
Family call	Tap Call Family	Dialer opens with the saved number; EasyCare does not call silently

Test	Steps	Expected result
Typed message	Type text and tap Send Message	SMS app opens with recipient and body filled in
Dictated message	Tap Speak Message and speak	Recognized text fills the message before SMS handoff
Quick emergency, location allowed	Allow location and send quick alert	SMS app opens; Maps link is included when a location fix is available
Quick emergency, location denied	Deny location and send quick alert	SMS app still opens with emergency text and no location; no crash
Spoken emergency	Dictate a custom emergency message	Detected text is used in the prepared SMS
Emergency call	Tap Call Emergency Services	Dialer opens with 112; call is not placed automatically
Medicine	Tap TAKEN	Item displays DONE ✓ and spoken feedback plays if enabled
Voice setting	Turn voice output off and visit screens	No Text-to-Speech feedback plays; controls still work
Text size	Select all three sizes	App text changes and important controls remain readable and usable
Missing external app	Test on a device without an SMS or dialer handler	EasyCare shows an error rather than crashing

14. Reset the app for another first-launch test

Resetting removes all saved EasyCare names, phone numbers, address, and settings.

On the device:

1. Open Android Settings.
2. Open Apps > EasyCare.

3. Open Storage.
4. Tap Clear storage or Clear data.
5. Open EasyCare again.

15. Troubleshooting

Voice recognition could not start

- Confirm microphone permission is enabled in Android Settings > Apps > EasyCare > Permissions.
- Check the internet connection.
- Confirm that voice recognition works in other apps on the phone.
- Wait for the current listening session to finish before tapping the microphone again.
- Continue with the large on-screen buttons if voice remains unavailable.

No location appears in an emergency message

- Turn on device Location.
- Allow precise or approximate location permission for EasyCare.
- Try outdoors or wait briefly for a location fix on a physical device.

An emergency SMS without a location is valid prototype behavior when no location is available.

The dialer or SMS app does not open

- Confirm a phone/dialer or SMS app is installed and enabled.
- Check that the saved phone number is not empty.

Setup appears again

Tap Save and Continue and allow the save to complete. Clearing app data or reinstalling the app intentionally returns EasyCare to first-launch setup.

Text or buttons do not fit on the screen

Ask the caregiver to select a smaller text size in EasyCare settings. Scroll vertically where the screen allows scrolling.